

β-actin兔多克隆抗体

β-actin Rabbit Polyclonal Antibody

本产品冰袋运输；-20℃保存，质保期为自收货后24个月内，-20℃保存时无需分装。

货号规格

货号	规格
LF202	100 μL
LF202S	20 μL



产品简介

本抗体为兔多克隆抗体,β-actin 的抗体是研究其在细胞内的分布和细胞骨架的有力工具。β-actin 是横纹肌纤维中的一种主要蛋白质，也是肌肉细丝及细胞骨架微丝的主要成分。它分布广泛，分子量为 42 kDa,具有高度保守性,在细胞中的表达相对稳定,因此常被用于校正系统的内参。


- 缓 冲 液10 mM PBS(pH 7.4), 50% glycerol
- 防 腐 剂0.02% sodium azide
- 应用种属人、大鼠、小鼠、猴、鸡、仓鼠、绵羊、非洲爪蟾
- 应用范围Western Blot(WB)、Immunohistochemistry(IHC)、Immunofluorescence(IF)
- 推荐浓度 (最适浓度需根据具体情况调试)
- Western Blot: 1:1,000 ~ 1:5,000;
- Immunohistochemistry: 1:1,000;
- Immunofluorescence: 1:1,000。

注意事项

1. 为了获得实验最佳效果, 请根据实验调整最佳工作稀释度;
2. 为了您的安全和健康, 请穿实验服并戴一次性手套操作;
3. 本产品仅限科研使用。



Immunohistochemical analysis of paraffin-embedded Heart tissue using Anti-β-Actin Rabbit pAb. Samples were incubated with primary antibody (1:1000) overnight at 4°C. A HRP-labeled anti-rabbit IgG (LF102) was used as the secondary antibody for 0.5 hour.



Immunofluorescent analysis of HeLa cells using Anti-β-actin Rabbit pAb at dilution of 1:1,000 and was visualized with Cy3 dye (Red).